

PocketChem™

UA
PU-4010

Every
Kidney
is worth
saving

National Kidney Foundation recommends:

- Measurement of UACR (Urine Albumin Creatinine Ratio) for detecting early stages of Chronic Kidney Disease (CKD*)
- Measuring UACR is as reliable as measuring the gold standard 24hr urinary albumin excretion
- Every Diabetic and Hypertensive should undergo an ACR test at least once in a year

Now you can measure UACR instantly and accurately in your clinic using SPOT urine sample

PocketChem™

UA PU-4010

Every Kidney is worth saving



US FDA
Approved

Made in
Japan

- Measures UACR with 95% correlation with Laboratory (Immunoturbidimetric) methods
ØRef: The Scandinavian Journal of Clinical and Lab investigation
- UACR is unaffected by variation in urine concentration / hydration status of the patient
- Test can be done any part of the day using SPOT urine sample
- Provides printed report for authenticity
- Also measures complete 10 routine urine parameters (Glucose, Protein, Bilirubin, Urobilinogen, pH, Specific Gravity, Blood, Ketones, Nitrites and Luekocytes)

ARKRAY Healthcare Pvt. Ltd.
701-702, Opulence, 6th Road, TPS III, Santacruz (E),
Mumbai - 400055, Maharashtra, India.

For more information contact:

022 3204 3204 | customersupport@arkray.co.in